

[Skip to content](#)[Curtir](#) {10} [Sitemap](#) [Contact us](#) [AGB](#) [Imprint](#)[Home](#) [Products](#) [About us](#) [Press & News](#) [Contact & Locations](#)You are here: [Home](#) » [Products](#) » [FREEWAY 014](#) » [Technical data](#)

[DIOR](#)  
[E-MAGIC PLUS](#)  
[APERTUS](#)  
[MAGICAL](#)  
[Genius MAGIC](#)  
[Genius TAXCOR I](#)  
[AMADEUS SUPERCROSS](#)  
**[FREEWAY 014](#)**  
[FREEWAY 035](#)  
[JOKER 014](#)  
[JOKER 035](#)

## FREEWAY™ 014 Paclitaxel Eluting Peripheral Dilatation Catheter

[Product information](#)   **[Technical data](#)**   [Order information](#)

### Technical data

FREEWAY™ 014 – latest, second-generation drug-releasing technology, 0.014" OTW infrapopliteal PTA balloon technology specifically designed for infrapopliteal interventions

#### FREEWAY™014: DEB

Design	Coaxial design - over the wire catheter (OTW)
Shaft diameter	dis. 3.1 F / medial 3.8 F / prox. 3.8 F
Balloon diameter	2.0 / 2.5 / 3.0 / 3.5 and 4.0 mm
Balloon length	40–150 mm
Usable catheter length (tip to strain relief)	120 cm and 150 cm
Guide wire diameter	0.014" (0.36 mm)
Coating	Hydrophilic coating on distal shaft
Balloon material	PA, Polyamid/Nylon
Balloon folding	Bi-fold for 2.0 mm Tri-fold for 2.5 to 4 mm
Balloon characteristic	Semi-compliant
Recommended introducer sheath	4 F
Nominal pressure	6 atm
Rated burst pressure	Diameter 2.0–2.5 mm: 16 atm/bar Diameter 3.0–4.5 mm: 14 atm/bar
Packaging unit	1 unit

[Go to top](#)[Home](#) [Products](#) [About us](#) [Press & News](#) [Contact & Locations](#)

Eurocor GmbH - In den Dauen 6a - 53117 Bonn - Germany **Phone:** +49 - 228 / 20 150 - 0 **Fax:** +49 - 228 / 20 150 - 5  
 © Copyright 2013. All rights reserved. Updated: 2017/03/31, Rev. No. 456671 (Version 20.0.0)  
 Eurocor GmbH is a wholly owned subsidiary of Opto Eurocor Healthcare Limited and is part of the Opto Circuits Group