

[Skip to content](#)Curtir [Sitemap](#) [Contact us](#) [AGB](#) [Imprint](#)[Home](#) [Products](#) [About us](#) [Press & News](#) [Contact & Locations](#)You are here: [Home](#) » [Products](#) » [FREEWAY 035](#) » [Technical data](#)

[DIOR](#)
[E-MAGIC PLUS](#)
[APERTUS](#)
[MAGICAL](#)
[Genius MAGIC](#)
[Genius TAXCOR I](#)
[AMADEUS SUPERCROSS](#)
[FREEWAY 014](#)
[FREEWAY 035](#)
[JOKER 014](#)
[JOKER 035](#)

FREEWAY™ 035 latest second-generation drug-releasing technology

[Product information](#)
 [Technical data](#)
 [Order information](#)

FREEWAY™ 035 DEB: Technical data	
Design	Bilumen design - over the wire catheter
Balloon diameter	4.0 / 5.0 / 6.0 / 7.0 and 8.0 mm
Balloon length	20–150 mm
Usable catheter length (tip to strain relief)	80 cm and 135 cm
Guide wire diameter	0.035" (0.91 mm)
Coating	Hydrophilic coating
Balloon material	PA, Polyamid/Nylon
Balloon folding	4-folding for 4 mm and 5 mm 5-folding for 6 to 8 mm
Balloon characteristic	Semi-compliant
Recommended introducer sheath	5 F for 4.0 to 6.0 mm Diameter 6 F for 6.0 and length ≥ 100 mm 6 F for 7.0 and 8.0 mm Diameter
Nominal pressure	6 atm
Rated burst pressure	Balloon length 20 / 40 / 60 mm Diameter 4–6 mm: 16 atm Balloon length 20 / 40 / 60 mm Diameter 7–8 mm: 14 atm Balloon length 80 / 100 / 120 / 150 mm Diameter 4 mm: 16 atm Balloon length 80 / 100 / 120 / 150 mm Diameter 5–6 mm: 14 atm Balloon length 80 / 100 / 120 / 150 mm Diameter 7–8 mm: 12 atm
Packaging unit	1 unit

[Go to top](#)

[Home](#) [Products](#) [About us](#) [Press & News](#) [Contact & Locations](#)

Eurocor GmbH - In den Dauen 6a - 53117 Bonn - Germany **Phone:** +49 - 228 / 20 150 - 0 **Fax:** +49 - 228 / 20 150 - 5
 © Copyright 2013. All rights reserved. Updated: 2017/03/31, Rev. No. 456671 (Version 20.0.0)
 Eurocor GmbH is a wholly owned subsidiary of Opto Eurocor Healthcare Limited and is part of the Opto Circuits Group