



Characteristics:

- Hydrophilic coating to ensure smooth introduction to navigate through tortuous anatomy
- Coil braided design provides outstanding pushability and kink resistance
- Radiopaque marker bands

Sneak Microcatheter Compatibility

Catheter	Size	Marker Bands	Total Length	ID	Embosoft size (µm)				
					100-300	300-500	500-700	700-900	900-1100
SNEAK 28	2,6F	1	130cm	0,71mm/0,028"	○	○	○	●	
SNEAK 24	2,4F	2	150cm	0,61mm/0,024"	○	○	○		
SNEAK 21	2,1F	2	150cm	0,53mm/0,021"	○	○	○		
SNEAK 17	1,9F	2	150cm	0,43mm/0,017"	○	○			